



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens High-Performance White Lab Coats

Product Code :
LBNY-00017-11000983



Description :

Wrinkle-resistant lab coat with five-button closure features 3 pockets for storage of tools and personal items.
Closure: five button
One-piece, lined, notched lapel collar
Pre-cure durable press finish
Care: light soil wash
Choose from regular or long length styles.
5 oz. performance blend poplin
80% polyester/20% combed cotton blend
Left breast pocket, barracked pencil stall, two lower pockets
Side vent openings
White.

Technical Specification :

Mens High-Performance White Lab Coats



Laboratory Instrument India

