



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mens High-Performance White Lab Coats

**Product Code :**  
LBNY-00017-11000983



**Description :**

Wrinkle-resistant lab coat with five-button closure features 3 pockets for storage of tools and personal items.  
Closure: five button  
One-piece, lined, notched lapel collar  
Pre-cure durable press finish  
Care: light soil wash  
Choose from regular or long length styles.  
5 oz. performance blend poplin  
80% polyester/20% combed cotton blend  
Left breast pocket, barracked pencil stall, two lower pockets  
Side vent openings  
White.

**Technical Specification :**

Mens High-Performance White Lab Coats



**Laboratory Instrument India**

