



Email : sales@laboratoryinstrumentindia.com

Product Name :
Top Dog 6 mil Die Cut 50 in. Apron

Product Code :
LBNY-00017-11000987



Description :

Designed to support zoned color-coding programs
One-piece construction limits the potential for foreign body contamination and reduces harborage points for pathogens
Resistant to common sanitation chemicals, and impervious to animal fats and degreasers.
Top Dog 6 mil Die Cut 50 in. Apron is reversible and designed to fit comfortably over clothes.
Reusable
FDA-compliant garments
Color-coded.

Technical Specification :

Top Dog 6 mil Die Cut 50 in. Apron



Laboratory Instrument India