



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Top Dog 6 mil Die Cut 50 in. Apron

**Product Code :**  
LBNY-00017-11000987



**Description :**

Designed to support zoned color-coding programs  
One-piece construction limits the potential for foreign body contamination and reduces harborage points for pathogens  
Resistant to common sanitation chemicals, and impervious to animal fats and degreasers.  
Top Dog 6 mil Die Cut 50 in. Apron is reversible and designed to fit comfortably over clothes.  
Reusable  
FDA-compliant garments  
Color-coded.

**Technical Specification :**

Top Dog 6 mil Die Cut 50 in. Apron



**Laboratory Instrument India**