



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Top Dog 6 mil X-Large Gown

**Product Code :**  
LBNY-00017-11000988



**Description :**

Combines the protection of sleeves and aprons into one garment  
Open back improves user comfort by reducing weight and providing ventilation  
Size x-large open-back gowns feature full-length cuffed sleeves and are available in a variety of colors.  
Color-coded.  
Reusable.

**Technical Specification :**

Top Dog 6 mil X-Large Gown



**Laboratory Instrument India**