



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Top Dog 8 mil Die Cut 50 in. Apron

**Product Code :**  
LBNY-00017-11000989



**Description :**

Ideal for workstations with sharp edges or jobs that require lifting, carrying, or the use of knives or scissors  
One-piece construction limits the potential for foreign body contamination and reduces harborage points for pathogens.

Designed to support zoned color-coding programs

Top Dog 8 mil Die Cut 50 in. Apron is reversible and designed to fit comfortably over clothes.

FDA-compliant garments

Thicker 8 mil material significantly increases resistance to punctures and tears

Reusable

Color-coded.

**Technical Specification :**

Top Dog 8 mil Die Cut 50 in. Apron



**Laboratory Instrument India**