



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Charcoal Lightweight Long Sleeve Flame Resistant T-shirt

Product Code :
LBNY-00017-11000991



Description :

Mens Lightweight Long Sleeve Flame Resistant T-shirt is a CAT 2 AR/FR, tee that features a heat seal triangle centered at the neck to show that the wearer has on an FR undergarment.

Fabric: Flame resistant, 7 oz. (235 g/m²) Interlock

Hazard: Flash Fire, Arc Flash

Wash Care: Home Wash

Blend: 100% Cotton.

Underarm gusset for improved breathability

Rib-knit crewneck collar

Single spade chest pocket.

Technical Specification :

Mens Charcoal Lightweight Long Sleeve Flame Resistant T-shirt



Laboratory Instrument India

