



Email : sales@laboratoryinstrumentindia.com

Product Name :
Phase Change Cooling Vest Charge Packs

Product Code :
LBNY-00017-11000992



Description :

Phase Change Cooling Vest Charge Packs are suitable for use with a variety of cooling vests.
Immediate and long-lasting cooling relief
Patented honeycomb ice pack design eliminates stiffness
Wash by hand with a mild detergent and warm water.
Packs maintain constant 64°F (18°C) for up to 4 hours
Quick recharge in a little as 10 minutes.

Technical Specification :

Phase Change Cooling Vest Charge Packs



Laboratory Instrument India