



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Phase Change Cooling Vest Charge Packs

**Product Code :**  
LBNY-00017-11000992



**Description :**

Phase Change Cooling Vest Charge Packs are suitable for use with a variety of cooling vests.  
Immediate and long-lasting cooling relief  
Patented honeycomb ice pack design eliminates stiffness  
Wash by hand with a mild detergent and warm water.  
Packs maintain constant 64°F (18°C) for up to 4 hours  
Quick recharge in a little as 10 minutes.

**Technical Specification :**

Phase Change Cooling Vest Charge Packs



**Laboratory Instrument India**