



Email : sales@laboratoryinstrumentindia.com

Product Name :
Type R Class 3 Long Sleeve T-Shirt

Product Code :
LBNY-00017-11000994



Description :

Full-sleeve coverage and maximum compliance from this breathable knit shirt with moisture wicking treatment.
For Use With (Application): Flagging Crew, Road Construction, Utility Crew, Survey Crew, Construction, Railway Workers, Warehouse/Distribution.

Color: Orange

No. of Pockets: 1

Size: 5X-Large.

Technical Specification :

Type R Class 3 Long Sleeve T-Shirt



Laboratory Instrument India