



Email : sales@laboratoryinstrumentindia.com

Product Name :
Type R Class 2 Black Front Performance Long Sleeve T-Shirt

Product Code :
LBNY-00017-11000995



Description :

Lightweight birds eye knit paired with Comfort Trim
For Use With (Application): Road Construction, Utility Crew, Survey Crew, Construction, Railway Workers, Baggage Handlers, Warehouse/Distribution.
Size: 5X-Large
Color: Orange

Technical Specification :

Type R Class 2 Black Front Performance Long Sleeve T-Shirt



Laboratory Instrument India