



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Type R Class 2 Two-Tone Mesh Vest with Reflective Binding

**Product Code :**  
LBNY-00017-11000996



**Description :**

Breathable mesh, plenty of storage and reflective binding accents.  
Material: Polyester.  
Color: Lime  
Size: 4X-Large to 5X-Large.

**Technical Specification :**

Type R Class 2 Two-Tone Mesh Vest with Reflective Binding



**Laboratory Instrument India**