



Email : sales@laboratoryinstrumentindia.com

Product Name :
PBI Bronze Aramid Fiber LTO Fire Protective Garment

Product Code :
LBNY-00017-11000999



Description :

Lightweight, one-piece FPG ensemble with snap-in removable liner
Elasticized waist for comfortable fit; articulating rapid-rescue strap
Double-stitched, lime 2 in.-wide reflective trim at chest, hip, wrist, and ankle
Reinforced elbows
Elasticized cuff with Velcro closure and over-the-thumb aramid fibers wristlets
Double-stitched pant cuffs with gray
Thermal liner of iced blue 7.5 oz. Semi-slick: 3.5 oz. filament/spun Nomex facecloth, and 4.0 oz. aramid fibers batt
Stand-up collar with closure
Outer shell of bronze 7.5 oz. aramid fibers ripstop fabric with enhanced water repellency
Front zipper with storm flap and Velcro closure
Zippered boot access panel on each leg.
Moisture barrier of 5 oz. with facecloth: 1.7 oz. Crosstech selectively permeable film laminated to 3.3 oz. pajama check woven substrate
Raglan sleeves for unrestricted movement, plus front zipper running from neck to crotch for quick donning/doffing
Two double-stitched 8 x 7.5 x 2.5 in. bellows pockets with reinforced bottoms and 5 in. flap
Double-stitched cushioned knees with 11 x 9 in. flap.

Technical Specification :

PBI Bronze Aramid Fiber LTO Fire Protective Garment



Laboratory Instrument India