



Email : sales@laboratoryinstrumentindia.com

Product Name :

Hot Rods Winter Liners: Shoulder Length

Product Code :

LBNY-00017-11001002



www.laboratoryinstrumentindia.com

Description :

Feature two built-in pockets to hold regular-size Hot Rods warming packs
Insulated ear barriers provide extra protection against wind, water and cold
Strip of reflective film above the back of neck offers visibility in low-light or dark hours of the day
Short elastic strip in the forehead area keeps liner snug to face
Two layers
Fleece lining
Flame-retardant until washed
100% cotton twill top
Color: Navy
Three layers
Fleece lining
Color: Green.
100% quilted nylon shell
2mm foam center layer.

Technical Specification :

Hot Rods Winter Liners: Shoulder Length

Laboratory Instrument India

