



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pullover Cooling Vest

Product Code :
LBNY-00017-11001003



www.laboratoryinstrumentindia.com

Description :

Super-absorbent, poncho-style workplace safety vest
Works without ice or freezing
Stays hydrated for days
100 percent cotton with super-absorbent cooling crystals
Expandable sides with adjustable buckles
Soak in cold water
One size fits all.

Technical Specification :

Pullover Cooling Vest



Laboratory Instrument India