



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pullover Cooling Vest

**Product Code :**  
LBNY-00017-11001003



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Super-absorbent, poncho-style workplace safety vest  
Works without ice or freezing  
Stays hydrated for days  
100 percent cotton with super-absorbent cooling crystals  
Expandable sides with adjustable buckles  
Soak in cold water  
One size fits all.

**Technical Specification :**

Pullover Cooling Vest



**Laboratory Instrument India**