



Email : sales@laboratoryinstrumentindia.com

Product Name :
Twill Action Back Coverall with Chest Pockets

Product Code :
LBNY-00017-11001004



Description :

This coverall is oversized for a roomy fit and can even be worn over clothes.
Finish: Wrinkle Resistant, Pre-cure durable press with Soil Release
Closure: Brass zipper, gripper at top and at lapel
Collar: One-piece, topstitched
Pocket: Two set-in front pockets, two patch hip pockets, rule pocket, and two breast pockets
Backside pleats offer additional room across the back of this coverall for easy movement and range of motion.
Vents also assist easy movement, so covering up won't limit you in getting your job done.
Action-back pleats for added mobility
Chest and hip pockets plus a handy rule pocket provide plenty of storage
Primary Closure: Zipper
Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash
Fabric: 7.25 oz. Twill
Blend: 65% Polyester/35% Combed Cotton
Part of the Zero Scratch Collection of worry-free uniforms to keep focus on work, not workwear
Side vent openings
Safety-stitched out seams and felled inseams
Four-needle band joins top and bottom
Gender: Male
Materials: Twill, Cotton/Polyester.
Product Collection: Zero Scratch
Fit Details: Fits true to size. Meant to be worn over clothes.
Regular length: 30 in. inseam for 5 ft. 7 in. to 5 ft. 11 in.
Long length: 32 in. inseam for 5 ft. 11 in. to 6 ft.3 in.

Technical Specification :

Twill Action Back Coverall with Chest Pockets



Laboratory Instrument India