



Email : sales@laboratoryinstrumentindia.com

Product Name :
Twill Action Back Coverall with Chest Pockets

Product Code :
LBNY-00017-11001004



Description :

This coverall is oversized for a roomy fit and can even be worn over clothes.

Finish: Wrinkle Resistant, Pre-cure durable press with Soil Release

Closure: Brass zipper, gripper at top and at lapel

Collar: One-piece, topstitched

Pocket: Two set-in front pockets, two patch hip pockets, rule pocket, and two breast pockets

Backside pleats offer additional room across the back of this coverall for easy movement and range of motion.

Vents also assist easy movement, so covering up won't limit you in getting your job done.

Action-back pleats for added mobility

Chest and hip pockets plus a handy rule pocket provide plenty of storage

Primary Closure: Zipper

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Fabric: 7.25 oz. Twill

Blend: 65% Polyester/35% Combed Cotton

Part of the Zero Scratch Collection of worry-free uniforms to keep focus on work, not workwear

Side vent openings

Safety-stitched out seams and felled inseams

Four-needle band joins top and bottom

Gender: Male

Materials: Twill, Cotton/Polyester.

Product Collection: Zero Scratch

Fit Details: Fits true to size. Meant to be worn over clothes.

Regular length: 30 in. inseam for 5 ft. 7 in. to 5 ft. 11 in.

Long length: 32 in. inseam for 5 ft. 11 in. to 6 ft.3 in.

Technical Specification :

Twill Action Back Coverall with Chest Pockets



Laboratory Instrument India