



Email : sales@laboratoryinstrumentindia.com

Product Name :
Jacket with Hood Snaps

Product Code :
LBNY-00017-11001005



www.laboratoryinstrumentindia.com

Description :

Provides Level C chemical protection
Features inner sleeve and cuffs, clear I.D. pocket, and reinforced elbows
Constructed of 0.40mm nylon/polyester lining coated with heavy-duty PVC
Detachable hood with heavy-duty draw string and cord lock snaps to jacket
Green.

Technical Specification :

Jacket with Hood Snaps



Laboratory Instrument India