



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Keystone Safety Wall Mountable Hair Net Dispenser

**Product Code :**

LBNY-00017-11001007



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Constructed from sturdy and durable High Impact PETG, making it 19 times stronger than regular acrylic.  
Reduce wastage of hairnets, bouffant and beard covers with this functional dispenser.  
Holds 100 pieces of hairnets/bouffant or 300 pieces of beard covers.  
Holds 100 hair nets/bouffant or 300 beard covers inside it and is wall mountable to better utilize your space.  
PETG Material  
Wall Mountable  
12 x 5 x 12 in.

**Technical Specification :**

Keystone Safety Wall Mountable Hair Net Dispenser



**Laboratory Instrument India**