



Email : sales@laboratoryinstrumentindia.com

Product Name :

Safety Laminated Polypropylene Shoe Covers - Aquasole

Product Code :

LBNY-00017-11001008



www.laboratoryinstrumentindia.com

Description :

Product Type: Controlled Environments Shoe Covers.

Material: Polypropylene

Color: White

Sole Type: Non-skid

Construction: Polypropylene

Size: X-Large

Weight (English): 12 lb.

Weight (Metric): 5.44 kg.

Technical Specification :

Safety Laminated Polypropylene Shoe Covers - Aquasole



Laboratory Instrument India