



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Dry Evaporative Cooling Vest

Product Code :
LBNY-00017-11001013



www.laboratoryinstrumentindia.com

Description :

Made of industrial grade quality nylon
Mesh side-panels provide additional ventilation and stretch
V-neck design with front zipper closure
Premium Dry Evaporative Cooling Vest offers immediate dry cooling relief that lasts up to three days.
Cooling material is treated with an anti-microbial that prevents mold growth
Dry technology keeps wearer completely dry and comfortable during the entire evaporative cooling process
Up to 3 days of cooling power to help maintain a healthy body temperature
Fill vest with 13 to 20 oz. (400 to 600 mL) of water
Reduces the effect of heat exhaustion and fights fatigue while working in the outdoor heat
Machine washable.

Technical Specification :

Premium Dry Evaporative Cooling Vest



Laboratory Instrument India

