



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Potentiometer

**Product Code :**  
LBNY-0003-1010063



**Description :**

Potentiometer

**Technical Specification :**

Constantan wire stretched along a metre scale, subdivided in cm and mm clamped to stout plated brass end plates.  
The terminating strips are provided with heavy duty 4 mm socket terminals and the whole is mounted on a polished hardwood baseboard.  
Comprising 0.56mm diameter constantan wire, scale 0-100cm and 100-0cm with a jockey and wooden base.



**Laboratory Instrument India**