



Email : sales@laboratoryinstrumentindia.com

Product Name :
Potentiometer

Product Code :
LBNY-0003-1010063



Description :

Potentiometer

Technical Specification :

Constantan wire stretched along a metre scale, subdivided in cm and mm clamped to stout plated brass end plates.
The terminating strips are provided with heavy duty 4 mm socket terminals and the whole is mounted on a polished hardwood baseboard.
Comprising 0.56mm diameter constantan wire, scale 0-100cm and 100-0cm with a jockey and wooden base.



Laboratory Instrument India