



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slotted Masses 10g

Product Code :
LBNY-0003-1010067



Description :

Slotted Masses 10g

Technical Specification :

Material: Made of chrome-plated steel, durable to use.

Accuracy: $\pm 2\%$ accuracy.

Convenient: Masses can be removed and replaced on the hanger.

Metric Weights: Set of metric slotted weights for weighing operations, for laboratory, teaching or home use.
Set of 10g.



Laboratory Instrument India