



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Ripple Tank

**Product Code :**  
LBNY-0003-1010068



**Description :**

Ripple Tank

**Technical Specification :**

To illustrate important wave motions like rectilinear and circular motion.  
This Ripple Tank is comprised of a molded frame with an acrylic window.  
Scratch resistant supported by adjustable legs.  
This dual purpose tank is a versatile apparatus allowing all types of waves and wave motion. With illuminator lamp and important accessories included, it forms a complete ready to use unit.  
The unique design and size allows it to be used on an overhead projector as well.  
Metal frame with an acrylic tray, four legs , a gently supplied for supporting a motor beam.  
Ripples are produced by a small 1.5V to 4.5V d.c motor mounted on a plastic beam of 240mm length, two spherical dippers, barriers (one spherical two long and one short)



**Laboratory Instrument India**