



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Fluid Pressure Apparatus Spouting Cylinder

**Product Code :**  
LBNY-0003-1010069



**Description :**

Fluid Pressure Apparatus Spouting Cylinder

**Technical Specification :**

To demonstrate that pressure increases with depth.

A solid metal cylinder 400 x 68mm diameter with three equal size orifices at different heights down one side.



**Laboratory Instrument India**