



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Press

Product Code :
LBNY-0003-1010070



Description :

Hydraulic Press

Technical Specification :

Adjustable hold-time at the bottom of the dead centre
Operation by foot pedal as well as by hand push buttons
Hydraulic ejector in the lower table
4 adapters for the piston head with pressure plates
Two pistons having capacity of 100 cubic centimetre and 20 cubic centimetre respectively and cross section area of 3 to 1 ratio.
Overall height is 240mm



Laboratory Instrument India