



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Stackable Plug Leads Connecting Wires

**Product Code :**  
LBNY-0003-1010072



**Description :**

Stackable Plug Leads Connecting Wires

**Technical Specification :**

Ideal for use where 4mm plug leads would be too bulky.  
The 2mm spring loaded plugs are connected by highly flexible PVC cable.  
4mm plug lead, nickel plated split hollow pins with maximum current 12A,  
Length: 150mm (red)



**Laboratory Instrument India**