



Email : sales@laboratoryinstrumentindia.com

Product Name :
Basic Metacentric Height Apparatus

Product Code :
LBNY-0005-1030002



Description :

To determine metacentric height of floating body and height variation with tilt angle. Tilt angle indicated by plumb bob on attached scale. Rectangular pontoon; centre of gravity can be moved side way by moving horizontal jockey weight, and vertically via adjustable vertical weight on mast.

Technical Specification :

Technical Detail:

Corresponding linear dimension: +90mm
Max. angle of heel: +13o
Pontoon dimensions: Length 350mm
Width: 200mm
Overall height: 475mm.



Laboratory Instrument India