



Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Phase Transformer Trainer

Product Code :
LBNY-0005-1050001



Description :

The trainer provides an experimental approach to the analysis of the behavior and characteristics of the single-phase transformer. Its operating features are deduced from typical industrial grade tests

Technical Specification :

2 dc outputs: 0-40V, 5A / 0-5V, 10A
2 ac outputs: 0-48V, 5A / 0-10V, 10A
2 ac/dc voltmeters: 10-50-100V
2 ac/dc ammeters: 1-5-10A
1 watt meter: 1-5-10A, 10-50-100V
Power supply: 220 V, 50/60 Hz
Resistive, capacitive and inductive loads.



Laboratory Instrument India