



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Lab

Product Code :
LBNY-0005-1060002



Description :

Digital Lab is used for error free transferring of different materials or data into digital formats via advanced network. This type of lab promotes real time communication between a group of people by eliminating possibility of time delay. Designed by efficient personnel, this lab consists of state of the art instruments and equipments for creating better working and learning environment. Special design of this laboratory helps to reduce its operation cost. We are a renowned manufacturer and supplier of excellent grade Digital Lab. Latest software is used to digitize video and audio files. Selected grade raw materials have been used to design this lab.

Features:

- Application of latest technical knowledge
- Precise diameter
- Simple to handle.

Technical Specification :

Digital Lab



Laboratory Instrument India

