



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Vessel

Product Code :
LBNY-0005-1070002



Description :

Once sterile, the UHT system is then placed in process mode and the into vessel fill mode. Sterilised product then fills the vessel. The UHT/filling system and the Sterile Vessel are first sterilised in separate processes. Time and temperature for the sterilisation process is set from the touchscreen.

The Sterile Vessel is controlled via a touchscreen mounted on the unit with modes for vessel filling, vessel emptying (bottling), sterilisation and cleaning. Vessel volumes range from 10 to 30 litres.

Features:

- Automated process operations
- Product agitation using a sterile gas supply
- 10-30-litre vessel volume
- Touchscreen display
- Steam sterilisable by external or integral supply
- Variable sterilisation protocol control from touchscreen display
- Continuous level sensing
- Clean-in-place capability
- Can be integrated into Armfield or other UHT/filling systems.

Technical Specification :

Sterile Vessel



Laboratory Instrument India