



Email : sales@laboratoryinstrumentindia.com

Product Name :
Car Diagnosis Box

Product Code :
LBNY-0002-640001



Description :

Car Diagnosis Box, Car Diagnosis Box Equipments, Car Diagnosis Box Tools, Car Diagnosis Box Tool Kits, Car Diagnosis Box Manufacturers, Car Diagnosis Box Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

Technical Specification :

Hardware is a little different, has 2 handles on the sides, and a fold prop for tool to stand, and some other minor hardware changes

Fault code read and clear

Current Data numerical and graphing

Actuation test

System Identification (ECU number)

Adaptation and coding dependent on vehicle and system

Printable Customer Relationship manager

Superior boot up and Connection Speed for Fast Service

Vehicle and Module auto detection(on most makes and models)

Built-in remote control and support availability via team viewer

Wireless internet connection for software updates and browsing the internet
Browse online for fault code/repair information by just clicking on the fault code
Wireless printing capability
8 inch IPS High Resolution Touch Screen
External memory storage compatibility
Car Diagnostic Scanner 8" IPS
LCD touch
Screen, 12-24 V
DC Input, USB2.0(1,4aster),Mini USB(Slave), RJ 45 Ethernet Cable
Connector, Interface Module Simulation



Laboratory Instrument India