




Email : sales@laboratoryinstrumentindia.com

Product Name :
Dimensional Metrology I Training Kit Three

Product Code :
LBNY-0005-1090002



Description :

The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. The test pieces used are stainless steel bushes with various inner and outer diameters. This dimensional metrology practice kit is designed for practical training in the metalworking trades. They were manufactured with the accuracy of CNC parts. All parts are clearly laid out and well protected on a plastic tray. The trays are stackable, providing for space-saving storage. The kit includes ten bushes, all of which differ to a minor degree in dimensions and each of which is individually marked.

Technical Specification :

Measurement exercises on a stainless steel bush
Instructional kit complete with test pieces and measuring aids
Practice kit for dimensional metrology in the metalworking trades
All parts clearly laid out and well protected on a plastic tray
10 Test pieces, each of different dimensions

Technical Data:

Test pieces: 6 lengths, 9 diameters measurable
Depth caliper gauge: 0...150mm
Three-point internal micrometer: d12...d16mm
Vernier caliper gauge: 0...200mm
Inside quick caliper: 10...30mm
Inside spring caliper: 125mm long
Pocket caliper gauge: 0...150mm

Internal micrometer: 25...50mm.



Laboratory Instrument India