



Email : sales@laboratoryinstrumentindia.com

Product Name :
Distillation Columns Unit

Product Code :
LBNY-0005-1100001



Description :

The equipment employs galvanically isolated intrinsically safe circuits and flameproof devices as appropriate, to enable safe distillation of standard test mixtures such as methylcyclohexane-toluene, methyl alcohol-water etc. The unit is a self-contained and fully instrumented distillation facility, suitable for practical laboratory work relevant to the teaching of unit operations.

The equipment consists of two interconnected units: a floor standing process unit and a bench-mounted control console.

Technical Specification :

Distillation Columns Unit



Laboratory Instrument India