



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Distillation Columns Unit

**Product Code :**  
LBNY-0005-1100001



**Description :**

The equipment employs galvanically isolated intrinsically safe circuits and flameproof devices as appropriate, to enable safe distillation of standard test mixtures such as methylcyclohexane-toluene, methyl alcohol-water etc. The unit is a self-contained and fully instrumented distillation facility, suitable for practical laboratory work relevant to the teaching of unit operations.

The equipment consists of two interconnected units: a floor standing process unit and a bench-mounted control console.

**Technical Specification :**

Distillation Columns Unit



**Laboratory Instrument India**