



Email : sales@laboratoryinstrumentindia.com

Product Name :
Modular Evaporator Family Unit

Product Code :
LBNY-0005-1100002



Description :

A laboratory evaporation system capable of being configured as rising or falling film, single or double effect. Integral USB interface for computer data logging and control. Contains vacuum pump, condenser and condensate vessel. Temperature controlled pre-heat stage. Adjustable re-circulation on each evaporation stage. Manual control console. Options to operate from external steam supply, or integral electrically heated pressurized hot water circulator.

Technical Specification :

Modular Evaporator Family Unit



Laboratory Instrument India