



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrical Boiler for Steam Production

Product Code :
LBNY-0005-1120002



Description :

The boiler water feeding is got from the mains by means of a embodied circulation pump rated at 0.55 kW and electrical valve. The automatic steam generator has a 9 kg/h production capacity, at 4.5 bar pressure. The most relevant protections and warning system are also considered, such as the level indicator, pressure switches, safety valves, manometer and warning lights. The heating of the water occurs by means of three resistors with a global electric power of 7 kW.

Technical Specification :

Technical Characteristics:

Working pressure: 4.5 bar
Steam production: 9 kg/h
General switch
Lights signaling the pump running and resistor supply
Water level indicator
Electric power: 7 kW, supplied by 3 resistors placed in the boiler
Manometer: 0 to 10 bar
Safety valve, pressure switch.

Laboratory Instrument India

