



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Automotive Scanner Tools and Car Diagnostic Scanner

**Product Code :**  
LBNY-0002-640002



**Description :**

Automotive Scanner Tools and Car Diagnostic Scanner

**Technical Specification :**

Equipped with a powerful for ideal management of workshop operations, including keeping all vehicle records, customer information, and other necessary images and data.

User interface through capacitive touch screen technology provides an easy & most intuitive use of a diagnostic tool and its multiple functions.

The high resolution TFT screen offers greater readability and convenient touch screen operation.

CPU : 2.1 GHz

LCD : ( G-scan-3) 10.1" TFT (1280 x 800 pixel)

Input Devices : Capacitive Touch Screen, 3 Hard keys (Power, Function 1, Function 2)

Automotive Scanner Tools and Car Diagnostic Scanner, Automotive Scanner Tools and Car Diagnostic Scanner Equipments, Automotive Scanner Tools and Car Diagnostic Scanner Tools, Automotive Scanner Tools and Car Diagnostic Scanner Tool Kits, Automotive Scanner Tools and Car Diagnostic Scanner Manufacturers, Automotive Scanner Tools and Car Diagnostic Scanner Suppliers from India, China, Kenya

---

# Laboratory Instrument India

