



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Two Way Hydraulic Valve

**Product Code :**  
LBNY-0005-1220002



**Description :**

Two hydraulic valve connected to the motorized pumping unit of the machine to activate a second frame. To activate alternatively two testing frames by using the same pumping unit. Installed on the pumping unit mod. Complete with protection case.

**Technical Specification :**

Two Way Hydraulic Valve



**Laboratory Instrument India**