



Email : sales@laboratoryinstrumentindia.com

Product Name :
Brinell Hardness Tester

Product Code :
LBNY-0002-640004



Description :

Brinell Hardness Tester

Technical Specification :

Brinell hardness tester that is ideal for wide range of Brinell loads from 62.5 kgf – 3000 kgf.

Designed with rugged construction to withstand harsh environments, the reliable Brinell hardness testing machine combines high rigidity and closed-loop load cell technology to ensure accurate and safe load applications.

Key features of the Brinell hardness tester include a clamping device for securing parts to the exchangeable anvil, an easy-to-use user interface to set up and operate the tester, a built-in keypad calculator that accurately presents HB hardness values, and an LCD screen.

Total Loads: from 10 to 3000kgf (From 98.07 to 29421N)

Brinell Hardness Tester, Brinell Hardness Tester Equipments, Brinell Hardness Tester Tools, Brinell Hardness Tester Tool Kits, Brinell Hardness Tester Manufacturers, Brinell Hardness Tester Suppliers from India, China, Kenya



Laboratory Instrument India

