



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turret Lathe

Product Code :
LBNY-0005-1230001



Description :

Features:

Enable repetitive production of same job requiring less involvement, effort and attention of the operator for Presetting of work speed and feed rate and length of travel of the cutting tools
Suitable for batch production or small lot production
Semi Automatic
More productive for quick engagement and overlapped functioning of the tools in addition to faster mounting and feeding of the job and rapid speed change
Possesses an axially movable index able turret in place of tail stock.

Technical Specification :

Turret Lathe



Laboratory Instrument India