



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Interactive PLC Trainer For Industrial Processes

**Product Code :**  
LBNY-0005-1260001



#### **Description :**

Palletizer: Palletize cases up to three layers  
Pick And Place: Place parts inside boxes through a three axes manipulator  
Automatic Warehouse: Transport, store and retrieve boxes from a rack  
Batching: Mix three primary colors (red, green and blue) in order to obtain a desired color  
Sorting: Transport cases from the entry bay to the elevators, sorting them by height.

#### **Features:**

Real time physics that closely emulates what happens in real life  
3D real time sound  
Friendly user interface  
Easy and fast to setup  
PLC Siemens of the Simatic  
3D realistic graphics  
Total interactivity with the systems.

#### **Technical Specification :**

Interactive PLC Trainer For Industrial Processes

---

# Laboratory Instrument India

