



Email : sales@laboratoryinstrumentindia.com

Product Name :
Interactive PLC Trainer For Industrial Processes

Product Code :
LBNY-0005-1260001



Description :

Palletizer: Palletize cases up to three layers
Pick And Place: Place parts inside boxes through a three axes manipulator
Automatic Warehouse: Transport, store and retrieve boxes from a rack
Batching: Mix three primary colors (red, green and blue) in order to obtain a desired color
Sorting: Transport cases from the entry bay to the elevators, sorting them by height.

Features:

Real time physics that closely emulates what happens in real life
3D real time sound
Friendly user interface
Easy and fast to setup
PLC Siemens of the Simatic
3D realistic graphics
Total interactivity with the systems.

Technical Specification :

Interactive PLC Trainer For Industrial Processes

Laboratory Instrument India

