



Email : sales@laboratoryinstrumentindia.com

Product Name :
Training System Level Control

Product Code :
LBNY-0005-1270001



Description :

Disturbance variables generated by electromagnetic proportional valve in tank outlet
Tank with overflow and graduated scale
Experimental unit for control engineering experiments
Level measurement by pressure sensor
Process schematic on front panel
Level control process with transparent tank
Speed-controlled pump.

Technical Specification :

Storage tank
Capacity: 3700mL
Level-controlled tank
Capacity: 1200mL

Pump

Power consumption: 18W
Max. head: 6m
Max. flow rate: 8L/min
Pressure sensor: 0...30mbar (0...300mm)
Process schematic with controller type selection
Proportional valve: Kvs: 0,7m/h.



Laboratory Instrument India