



Email : sales@laboratoryinstrumentindia.com

Product Name :
Training System Pressure Control

Product Code :
LBNY-0005-1270002



Description :

Electronic pressure sensor
Solenoid valve to generate disturbance variables
Pressure control in a tank
Software-based controlled system simulation
Experimental unit for control engineering experiments
Process schematic on front panel
Speed controlled diaphragm gas pump.

Technical Specification :

Diaphragm gas pump
Max. flow rate: 3L/min
Max. negative pressure: 250mbar abs.
Max. positive pressure: 1bar
Pressure tank
Operating pressure: 1bar
Capacity: 400mL
Pressure control range: 0...1bar
Solenoid valve: Kvs: 0,11m³/h
Max. pressure: 10bar
Pressure transducer: 0...1bar
Manometer: 0...1bar.



Laboratory Instrument India