



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Control Valve Characteristics Apparatus

**Product Code :**  
LBNY-0005-1280002



**Description :**

The Rotameter is used for the flow measurement. The pressure at the valve inlet is measured with the help of manometer. The air regulator and pressure gauge is provided for regulating air supply. The setup is specially designed for various experimental and demonstrative capabilities. The Stainless Steel water tanks with pump are provided for continuous water circulation. These units along with necessary piping and fitting are assembled in a painted MS Structure. The product consists of Pneumatic control valves of Linear and Equal % type.

**Technical Specification :**

**Experimentation:**

Flow Co-efficient (Cv) of Control Valve.  
Hysteresis of Control valve.  
Installed Characteristics of Control Valve.  
Rangeability of Control valve.  
Inherent Characteristics of Control Valve.



**Laboratory Instrument India**

