




Email : sales@laboratoryinstrumentindia.com

Product Name :
Vertical Bench Top Frame

Product Code :
LBNY-0005-1290002



Description :

Vertical Bench Top Frame A set of hangers and weights are supplied which enable the systems to be loaded. A variety of weights are supplied which gives accurate loading. The kit comprises a set of single, double and triple pulleys that allow a simple or compound pulley system to be constructed. Hooks on each pulley set allow the easy interconnection of pulleys to generate different system. A Vertical Bench Top Frame has a backboard attached with pegboard holes throughout. Vertical Bench Top Frame allows easy and quick installation of the component parts from the kit supplied. Cords are supplied that interconnect the pulleys, but also to use to run over each pulley arranged. All component parts are supplied neatly and safely into a carry case.

Technical Specification :

Various cord lengths
Hangers and Weights provided
Single, double and triple pulleys set
Vertical frame and backboard
Tape measure
Size of unit 400mm x 650mm x 400mm.

Experimental Capabilities:

Mechanical Advantage
Efficiencies
Pulley Ratio.



Laboratory Instrument India