



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Nozzle Pressure Distribution Module

**Product Code :**  
LBNY-0005-1320001



**Description :**

All three nozzles have axial static pressure tapings allowing the approach, throat and divergent section pressures to be measured. The variation of pressure ratio and mass flow may be investigated for all three nozzles.

Two convergent-divergent nozzles of with the same throat diameter but different discharge area and a single convergent nozzle having the same diameter are supplied and fit in the common test section.

**Technical Specification :**

Nozzle Pressure Distribution Module



**Laboratory Instrument India**